



All products are CLIA Waived and either 510(k) Cleared or Exempt.

hGC Pregnancy Test

Rx = For Prescription Use Only
 OTC = For Over-the-Counter Use

For the qualitative detection of human chorionic gonadotropin in urine or serum to aid in the early detection of pregnancy.



- Results within 5 minutes
- **99% Accurate
- **Detects pregnancy 6 days sooner than the date of the expected period

SKU*	Units	Cutoff	Format	Rx/OTC
GAHCG-101a-50, 100	50/100 tests/kit	25 mIU/mL	Strip	Rx
GAHCG-101a-H5, 10	5/10 tests/kit			OTC
GAHCG-102a	25 tests/kit	25 mIU/mL	Cassette	Rx
GAHCG-102a-H1,2	1/2 tests/kit			OTC
GAHCG-103a-H1,2	1/2 tests/kit		Midstream	OTC
GAHCG-103m**-H1,2	1/2 tests/kit	10 mIU/mL		OTC

*For different test sizes, please inquire with a sales representative

LH Ovulation Rapid Test

For the qualitative detection of luteinizing hormone (LH) in urine to aid in the detection of ovulation.



- 99% Accurate
- Positive results indicate ovulation should begin within the next 24 - 36 hours
- Results in 10 minutes



Digital Ovulation (LH) Monitor

SKU	Units	Cutoff	Format
GALH-101a	100 Tests/Kit	25 mIU/mL	Strip
GALH-101d		30 mIU/mL	
GALH-101b		40 mIU/mL	
GALH-102a	25 Tests/Kit	25 mIU/mL	Cassette
GALH-102d		30 mIU/mL	
GALH-102b		40 mIU/mL	
GALH-103a	1/5 Test(s)/Kit	25 mIU/mL	Midstream
GALH-103d		30 mIU/mL	
GALH-103b		40 mIU/mL	

FSH Menopause Test

For the qualitative detection of Follicle Stimulating Hormone (FSH) level in urine to evaluate the onset of menopause in women.

- 99% Accurate
- Requires Serial Testing; 2 tests 1 week apart
- Results in 5 minutes




SKU	Units	Cutoff	Format
GAFSH-101a	100 Tests/Kit	25 mIU/mL	Strip
GAFSH-102a	25 Tests/ Kit		Cassette
GAFSH-103a	1 Test/ Kit		Midstream

UTI Test Strip

For the qualitative detection of leukocytes (white blood cells) and nitrite in urine as an aid in the screening of a Urinary Tract Infection (UTI).

- Results in 2 minutes





SKU	Units
UTI-2T-H1	1 Test
UTI-2T-H3	3 Tests
UTI-2T-H6	6 Tests

Vaginal pH Rapid Test

For the qualitative detection of pH in female vaginal swab specimens to aid in the diagnosis of bacterial vaginosis.

- pH Range 3.8-7.0
- Results under 2 minutes



SKU	Units	Format	Rx/OTC
VPH-501a	25 Tests/Kit	Strip	Rx
VPH-501a-H2,5	2/5 Tests/Kit		OTC
VPH-502a	25 Cassettes/Kit	Cassette	Rx
VPH-502a-H1,2	1/2 Cassettes/Kit		OTC